



Buffer\$ 2.0 README File

Running Buffer\$ from the CD

Buffer\$ can be run from the CD by clicking on the file *Buffer\$ 2.0.xls*. When opening the file, you should see a screen like this:




Select the Enable Macros button. The screen will then look like this:



We recommend viewing the Instructional Videos to quickly learn how to use the tool and its various components.

Running Buffer\$ from a hard drive.

You may decide to run Buffer\$ from your hard drive to have the program run quicker and to avoid having the CD in CD drive. Create a folder called Buffer\$ on your hard drive and copy all of the CD files to this folder.



The program can be run by selecting the Buffer\$ 2.0.xls file on your hard drive. Since the videos were also copied, they will be accessible in the Buffer\$ program.

If you would like a shortcut on your desktop for opening Buffer\$, drag the Buffer\$ 2.0 Shortcut on to the desktop. Once it is on the desktop, you can double-click on the shortcut to start the program.

